

Esmarch Bandages



Esmarch Bandages

Length: 5m, autoclavable up to 121°C

Width: 6cm	Material: rubber, not latex free	REF. 20-20-060
Width: 8cm	Material: EPDM, latex free	REF. 20-20-080
Width: 10cm	Material: EPDM, latex free	REF. 20-20-100
Width: 12cm	Material: EPDM, latex free	REF. 20-20-120
Width: 15cm	Material: EPDM, latex free	REF. 20-20-150



Esmarch Bandages made of Silicone

Length: 4m, autoclavable up to 134°C

Width: 8cm	Material: Silicone, latex free	REF. 20-22-080
Width: 10cm	Material: Silicone, latex free	REF. 20-22-100