

bluestar

Face mask devised for Non invasive assisted and/or controlled Ventilation (CPAP and Bilevel). The cushion guarantees an optimal air-leak control and high comfort level. Due to the special contour adaptation of the cushion, skin break-down is minimal. The mask also comes with a special connector with ports for O₂, EtCO₂ and/or airways pressure monitoring. Single patient use.



SIZE INFORMATION

Code	Description	Size	Weight	Dead space	Connector Ø
KPM-001-L	NIV High Pressure Face mask	Ultra Large	91,0 gr	290 ml	External 27 mm Internal 22 mm
KPM-001-M	NIV High Pressure Face mask	Large	79,5 gr	240 ml	External 27 mm Internal 22 mm
KPM-002-M	NIV High Pressure Face mask	Medium	70,0 gr	220 ml	External 27 mm Internal 22 mm
KPM-003-M	NIV High Pressure Face mask	Small	65,0 gr	175 ml	External 27 mm Internal 22 mm
KPM-004-M	NIV High Pressure Face mask	Child	46,0 gr	85 ml	External 27 mm Internal 22 mm
NIV	Face mask connector		19,0 gr	17 ml	External 27 mm Internal 22 mm