


# Pediatric Breathing Systems

## Anesthetic-System

autoclavable

- consisting of:

- 2 x 10mm Silicone Respiration Tubing, 110cm long, 22mm/10mm Cuff
- 1 x 10mm Silicone Respiration Tubing, 110cm long, 22mm/22mm Cuff
- 1 x Y-Connector made of TPX
- 1 x Rebreathing Bag 500ml, made of Silicone




Ref. no. 62-11-200

ditto, but with Rebreathing Bag made of Latex

Ref. no. 62-11-100

- special versions upon request

## Silicone Catheter Mount Tubings




Catheter Mount Tubing 10mm I.D. / 15mm I.D., 10,5cm long, autoclavable

Ref. no. 65-00-105

Catheter Mount Tubing 10mm I.D. / 10mm I.D., 9cm long, autoclavable


Ref. no. 65-00-090

## Connecting Pieces



Y-Connector, swivelling, with two Silicone Catheter Mount Tubings for Children with 15mm Connectors

Ref. no. 65-20-100



Y-Connector, swivelling, with Suction Port, Silicone Suction Cap and two Silicone Catheter Mount Tubings for Children with 15mm Connectors

Ref. no. 65-20-200

**VBM** MEDIZINTECHNIK GmbH