

PARKS

Pocket Doppler *with built-in speaker*



Model 841-A

This Doppler is successor to the Model 841, our most popular pocket Doppler. It is an updated version which allows plug-in probes so you can use different types of probes or just for convenience in replacing defective probes. It also features a timer which will shut off the unit in five minutes. Pushing the knob in turns the Doppler on or off. The built-in speaker is convenient for teaching or to let the patient hear.

We recommend this Doppler be used where portability is of prime importance. It clips to your belt or fits in your coat pocket. While it was made for vascular surgeons, many hospitals do use it on the wards. Small instruments, however, have a tendency to disappear.

If you need a Doppler with you for making rounds with residents and want to listen to the posterior tibial, common femoral, digital or supraorbital arteries, this little Doppler will put on an impressive performance. The nominal 8 MHz. frequency will give a clean, clear sound when arterial flow is diminished or when you are checking venous flow. It is slightly less sensitive than our full size Dopplers like Model 811-B. The limitation is due to available battery power.

We have provided an auxiliary output on the Model 841-A so you can connect a stethoscope-type earphone if you wish.

Plugging in the earphone shuts off the loudspeaker automatically. A stethoscope earphone is available at a small extra cost.