

VBM ROLL-ON CUFFS and ACCESSORIES

Silicone Roll-On Cuffs



Size	Colour Code	Extremity Circumference	Ref.No.
0	grey	16-21 cm	21-95-710
1	white	20-28 cm	21-95-711
2	light-blue	26-33 cm	21-95-712
3	yellow	31-40 cm	21-95-719
4	red	38-49 cm	21-95-722
5	green	47-60 cm	21-95-727
6	dark-blue	58-70 cm	21-95-729
7	transparent	67-80 cm	21-95-730

Accessories for Roll-On Cuffs



Manometer with 1m Connection Tubing, Ref.No. 21-95-300



Tape-Measure, colour coded, Ref.No. 21-95-111



Brake-Block for Roll-On Cuffs, Ref.No. 21-95-200

Roll-On Cuff Set



Roll-On Cuff Set, consisting of Roll-On Cuffs size 1-4, Manometer with 1m Connection Tubing, Tape-Measure, two Brake-Blocks and transparent Case, Ref.No. 21-95-400