

VBM SINGLE CUFFS

Single Cuffs made of fabric with silicone bladder, completely washable



Ref. No.	Description	Cuff size	
		cm	inches
		length x width	length x width
20-44-700	Babies	30 x 3.0	12.0 x 1.0
20-44-710	Children	40 x 6.0	16.0 x 2.5
20-44-711	Adults, arm	60 x 7.5	23.5 x 3.0
20-44-712	Adults, arm	70 x 11.0	27.5 x 4.5
20-44-719	Adults, arm	100 x 11.0	39.5 x 4.5
20-44-721	Adults, thigh, <i>conical</i>	70 x 11.0	27.5 x 4.5
20-44-722	Adults, thigh, <i>conical</i>	100 x 11.0	39.5 x 4.5
20-44-727	Adults, thigh, <i>conical</i>	120 x 11.0	47.0 x 4.5
20-44-729	Adults, thigh, <i>conical</i>	135 x 11.0	53.0 x 4.5
<i>LOW PRESSURE CUFFS, extra wide</i>			
20-44-819	Adults, arm / thigh	110 x 14.5	43.5 x 6.0
20-44-827	Adults, thigh, <i>conical</i>	130 x 14.5	51.0 x 6.0
20-44-829	Adults, thigh, <i>conical</i>	145 x 14.5	57.0 x 6.0

Silicone Single Cuffs with blue strap completely washable and autoclavable up to 134°C (273°F)



Ref. No.	Description	Cuff size	
		cm	inches
		length x width	length x width
20-94-710	Children	44 x 6.0	17.5 x 2.5
20-94-711	Adults, arm	55 x 9.0	21.5 x 3.5
20-94-712	Adults, arm	66 x 9.0	26.0 x 3.5
20-94-722	Adults, thigh	90 x 10.5	35.5 x 4.0
<i>LOW PRESSURE CUFFS, extra wide</i>			
20-94-819	Adults, arm / thigh	92 x 16.0	36.0 x 6.5
20-94-825	Adults, thigh	92 x 16.0	36.0 x 6.5

Silicone Single Cuffs with white velcro fastening completely washable and autoclavable up to 134°C (273°F)



Ref. No.	Description	Cuff size	
		cm	inches
		length x width	length x width
20-84-700	Babies	24 x 5.0	9.5 x 2.0
20-84-710	Children	35 x 7.0	14.0 x 3.0
20-84-712	Adults, arm	75 x 12.0	29.5 x 5.0
20-84-722	Adults, thigh	100 x 12.0	39.5 x 5.0
20-84-727	Adults, thigh	115 x 15.0	45.5 x 6.0
20-84-729	Adults, thigh	125 x 15.0	49.5 x 6.0

Standard Cuff

VBM

CE