

740

Monitor from CAS

Created with emergency medical responders in mind, this multi-parameter monitor sets a new standard for accuracy, reliability and ease of use.

The CAS 740 monitor offers leading-edge technology in a rugged, lightweight, portable unit.

CAS Medical's breakthrough MAXNIBP™ technology reduces the problems of motion in blood pressure measurement. The 740 also lets you monitor oximetry and temperature, either continuously or on a "spot-check" basis, eliminating the need for additional equipment. Oximetry choices include Masimo SET® pulse oximetry, which is designed to be accurate during motion and low perfusion.

Yet the 740 weighs only 3 pounds; operates on battery, AC line or emergency vehicle power; and is simple to program. It can be fixed to a swivel wall mount designed specifically for emergency vehicles or carried in a protective, durable case.



Swivel wall mount



Carrying case

EMS Applications



 **CAS MEDICAL SYSTEMS, INC.**
TECHNOLOGY APPLIED TO MEDICINE

CAS Medical Systems, Inc. Telephone: 800 227 4414 or 203 488 6056. www.casmed.com