


# ANESTHESIA

## *crystal*

Respiratory Anesthesia face mask with inflatable air cushion.  
High transparency allows good visibility.  
Available in different scents. Single use.



### SIZE INFORMATION

Code	Description	Size	Weight	Dead space
<b>KM-200</b>	Respiratory/Anesthesia mask	Ultra Large	59,5 gr	310 ml
<b>KM-201</b>	Respiratory/Anesthesia mask	Large	55,5 gr	250 ml
<b>KM-202</b>	Respiratory/Anesthesia mask	Medium	51,0 gr	220 ml
<b>KM-203</b>	Respiratory/Anesthesia mask	Small	37,5 gr	180 ml
<b>KM-204</b>	Respiratory/Anesthesia mask	Child	27,0 gr	85 ml
<b>KM-205</b>	Respiratory/Anesthesia mask	Infant	17,5 gr	45 ml
<b>KM-206</b>	Respiratory/Anesthesia mask	Neo-Nate	14,5 gr	25 ml
<b>KM-207</b>	Respiratory/Anesthesia mask	Pre Neo-Nate	11,0 gr	20 ml