


headgear



Soft cap for face and nasal masks for
Non Invasive Ventilation

SIZE INFORMATION

Code	Description	Size
KPM-P001	Headgear for Face masks	Large
KPM-P002	Headgear for Face masks	Small
KM-P003	Headgear for Nasal masks	Large
KM-P004	Headgear for Nasal masks	Small

tri-gym

System with three coloured balls for flow breathing exercises IN (inspiration) and OUT (expiration). During the breathing exercises, it is also possible to have PEP or PIP resistances, by modifying the special IN/OUT valves.



SIZE INFORMATION

Code	Size
KM-805	Standard